Promedon
People + Innovation

PROMEDON
Av. Gral. Manuel Savio s/n
Lote 3 • Manzana 3 • (X5925XAD)
Pque Industrial Ferreyra, Cba • Argentina

www.promedon.com

CE 0197

M-20-15-0114 (05)

n&s

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
USER MANUAL

DESCRIPTION

The N&S Testicular Prosthesis has been designed to replace, through surgery, the natural testicle in shape and size. It is a solid implant made of elastomer of low-hardness silicone, which is very similar to the natural testicle.

PRESENTATION

IMPLANTS		
MODEL	A (mm)	B (mm)
T-Small	31	23
T-Medium	37	28
T-Large	42	32
T-XLarge	47	34

Sterilization Method: Ethylene Oxide.

INDICATIONS

Designed as a cosmetic replacement, N&S Testicular Prosthesis is advisable for those patients who had suffered from testicular agenesis, or when the testicle has been removed through surgery due to several pathologies.

CONTRAINDICATIONS

The prosthesis implant is not recommended if there is infection in some part of the body, particularly if it is genital or in the urinary tract.

N&S Testicular Prosthesis cannot be used in patients with neoplasia or transscrotal infection.

WARNINGS

Patients who undergo surgical operations are liable to complications during and after the surgery. Surgeries associated with the use of testicular implants entail risks or potential complications. Thus, prior

to surgery, physicians should inform patients about possible complications related to the use of this prosthesis. This product has been designed for single use.

Testicular implants in children is only recommended for cosmetics purposes.

Patients must be taught to distinguish the prosthesis from the natural testicle through self-examination.

Removal of the prosthesis is advisable in the case of surgical, physical or psychological problems.

PATIENT INFORMATION

It is the surgeon's responsibility to inform the patient and/or his representatives, before surgery, about all possible complications related prosthesis implantation.

Promedon and its distributors delegate to the surgeon the responsibility of informing the patient about all possible risks and advantages related to prosthesis implantation and usage.

PRECAUTIONS

Careful patient selection is essential, as well as a thorough diagnostic study before the surgery.

The prosthesis must be checked prior to and during the surgery, so as to monitor the structural integrity of the implant.

The prosthesis should be handled carefully, avoiding pointed, toothed, or sharp instruments, as any nick or damage may be the cause for subsequent complications of the implant.

Dirt, fingerprints, talc or any other substances that can contaminate the surface of the prosthesis may be the cause of reactions to foreign objects. Extreme preventive measures must be taken to avoid contamination.

Any hole or mark in the prosthesis is a possible cause of failure, since it may serve as a surface to keep materials that may cause reactions to foreign objects or infections in the patient.

The prosthesis is supplied sterile and pyrogen free in a double-wrap packaging system.

This packaging is contained into a cardboard box wrapped with a thermocontribute foil which must be undamaged. IF IT WAS DAMAGED WHEN HANDLING, THE PROSTHESIS SHOULD NOT BE IMPLANTED.

Operating room conditions: In general, operating rooms meet the requirements for testicular implant surgery. Nevertheless, we strongly recommend checking the following parameters:

- Operating room asepsis.
- Adequate training for the surgical support personnel.

COMPLICATIONS

Possible complications associated with the use of this prosthesis must be discussed with the patient prior

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to surgery.

Complications which may result from the use of this prosthesis include the risks associated with the medication and methods utilized in the surgical procedure, as well as the patient's degree of intolerance to any foreign body. Some complications may demand prosthesis removal.

Pain or fever indicating infection may appear after the implantation.

Infections that do not respond to antibiotic therapy demand prosthesis removal.

Dermic necrosis or wound opening may appear as a result of: an inadequate tension of the skin covering the implant, a trauma of skin surface during surgery or an inadequate circulating inhibition of the tissue. Subsequent extrusion of the implant may be necessary after this.

Post-surgical hematoma, seen as a swelling and tissue color fading, may cause extrusion of the device, if it not treated appropriately.

Excessive fluid accumulation around the implant may be produced after surgery, as a result of traumas.

Prosthesis size, wrong placing or migration may cause unsatisfactory cosmetic results.

The post-surgical formation of a fibrous tissue capsule surrounding the testicular prosthesis is a normal physiological reaction to the implantation of a foreign body.

If the patient feels discomfort, pain, throbbing of heart or prosthesis movement, the implant must be removed. Surgeons should inform **Promedon S.A.** or **Promedon's** distributor about any complication associated with the use of this prosthesis.

SURGICAL APPROACH

Although the recommended approach is via inguinal, the surgeon should choose the most convenient approach for each particular case.

In the case of patient with unilateral agenesis, the date of implantation should be defined by the surgeon. When orchietomy is required, ablation and implant are performed during the same surgery.

In those patients with testicular tumours, cord torsion, atrophic testicles or elderly patients with inguinoscrotal hernia, just one surgical act is performed.

In the case of patients with traumatism, it is advisable that the surgeon try a reconstruction without implant. However, if the results is negative, the surgeon should implant the prosthesis.

In case of castration for prostatic carcinoma, the surgeon should choose the most convenient approach for each patient, but it is common to implant the prosthesis during the same surgery.

It is for the own judgement of the surgeon to decide to fix the prosthesis or not.

In case he fixes it, he has to pass three suturing material of his preference through whichever point of the prosthesis, since its solid design avoids leakage risks. However, we suggest that a minimal distance of 5 mm from the extreme of the prosthesis should be respected so as to avoid the prosthesis breakage.

The bibliography mentioned bellow could help surgeons to get acquainted with the techniques.

REFERENCES

- Michael Beer, M.D., and Robert Kay, M.D., Testicular Protheses, Urologic Clinics of North America - Vol. 16, N° 1, February 1989, pages 133 to 138.
- Joseph Ortengerg, M.D., and Robert G. Kujpper, M.D., "The Why and How of Synthetic Replacement Testicles", Contemporary Urology, October 1991, pages 23 to 32.

SYMBOLS USED IN THE LABELS

	CATALOGUE NUMBER
	BATCH NUMBER
	SEE INSTRUCTIONS FOR USE
	EXPIRATION DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE: STERILIZATION METHOD ETHYLENE OXIDE
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

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