|                                |     |
|--------------------------------|-----|
| 1. French - Français.....      | 4   |
| 2. English.....                | 23  |
| 3. Danish - Dansk.....         | 42  |
| 4. Dutch - Nederlands.....     | 61  |
| 5. Finnish - Suomi.....        | 80  |
| 6. German - Deutsch.....       | 99  |
| 7. Italian - Italiano.....     | 118 |
| 8. Portuguese - Português..... | 137 |
| 9. Spanish - Español.....      | 156 |
| 10. Swedish - Svenska.....     | 175 |

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
# tensi+

## DEVICE DESCRIPTION

Tensi+ was designed to stimulate the posterior tibial nerve.

**1. Positive electrode**

**2. (+) button**, to increase the stimulation intensity

**3. Micro-USB charging port and charge indicator**

**4. LEDs** showing the stimulation intensity

**5. Retention strap loops**

**6. (-) button**, to reduce the stimulation intensity

**7. Battery charge level indicator**

**8. On/Off button** : press and hold for 1 sec to switch the device on or off

**9. Negative electrode**

**Precautions for use**

Tensi+ must only be used with the original accessories (retention strap, silicone electrodes). The device must always be placed along the posterior tibial nerve on the inner side of the (right or left) ankle.

The device must never be placed in another area. Place the device in the correct direction, with the on/off button at the bottom.

**Warnings**


Do not use Tensi+ for other indications.

Simultaneous use with a high-frequency surgical device may cause burns under the Tensi+ electrodes.

Tensi+ is a medical electronic device requiring special precautions for electromechanical compatibility. This device should not be used in an environment with strong electromagnetic interference. Using Tensi+ near (around 1 m/3 ft. from) a microwave, short-wave or ultra-short-wave device (microwave oven, Wi-Fi terminals or modules, mobile phones, Bluetooth systems, etc.) may alter the stimulation settings.

**Undesirable effects**

If you experience pain or discomfort, stop using the product and contact your doctor. Use maximum intensity treatments with caution. Do not exceed your comfort level.

### 3. ABOUT TENS TECHNOLOGY

Transcutaneous Electrical Nerve Stimulation, more commonly known as TENS, is a non-medicinal, non-invasive technique. It uses electrodes on the skin to deliver electrical pulses close to an innervated area to affect nerve transmission. Tensi+ targets the posterior tibial nerve, part of which originates close to the nerve fibres controlling the bladder. Stimulating this nerve modulates messages sent between the bladder and the brain, improving the central nervous control of urination. Clinical studies have shown that one daily 20-minute treatment reduces\* the symptoms of overactive bladder (urinary urgency, urge incontinence, nocturia or pollakiuria) and improves quality of life.

*\*A decrease in symptoms was observed from 1 month of treatment.  
It is recommended to continue treatment for at least 2 months*

### 4. TECHNICAL DESCRIPTION

Tensi+ is a Transcutaneous Electrical Nerve Stimulator (TENS) placed directly on the posterior tibial nerve to deliver pulses of an adjustable intensity. The Tensi+ device has a lifespan of 2 years according to IEC 60601-1-11.

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|                |                       |
|----------------|-----------------------|
| Output current | 0 to 50 mA (+/- 10%)  |
| Frequency      | 10 Hz (+/- 20%)       |
| Pulse width    | 200 $\mu$ S (+/- 20%) |
| Rated current  | 15 mA                 |

The values are valid under a resistance of 1000 Ohm +/- 10%.

**b. Assembling the device**

The device is already assembled with 2 textile ankle straps. To adjust the device to your ankle, pull on either the large or small strap.

If the pre-assembled device is too loose, remove the small strap and only use the large strap. The strap keeps the device comfortably in place for use.

Depending on your preference, position the plastic strap buckle on either the inner side or outer side of your ankle. 2 pre-assembled electrodes are provided with the device.

Make sure they are installed correctly (see section 8.a) before using the device.

**Assembly with 2 straps****Assembly with 1 strap**

### c. Positioning the device

#### **Step 1:** Clean your skin

Always use the device on clean skin. Any oil, cream, lotion, dust or other substance on the skin could affect the electrical pulses.


Clean the area to be treated (inner side of your ankle and Achilles tendon) with a damp cloth. Using the device on slightly damp skin improves the electrical pulses.

#### **Step 2:** Add gel to the electrodes (recommended)

Place a small drop of gel on each electrode, just before applying Tensi+ to the site indicated in the next section.

#### **Step 3:** Place the device on the ankle

Install the device with the power off on either your right or left ankle, as shown below.


1. The lower electrode is placed just below and behind the medial malleolus (inner side of the ankle).

2. The device is then placed in the same direction as the leg.

Once the device is in place, pass the free end of the strap through the plastic buckle and tighten to adjust to your ankle.

For your comfort, do not pull too tight.

**Please note:**

The device will tell you if you exceed 10 mA by emitting 2 short beeps. This is simply for information purposes.

To ensure that the device is positioned correctly on the tibial nerve, you can increase the intensity until your big toe contracts.

If the device is detached from the skin for more than one second, the stimulation level immediately drops to zero to keep the user safe and all 6 LEDs flash. You will therefore need to adjust the device again.

If you feel tingling under any of the electrodes, it means that the electrode is not properly in contact with your skin. Check that it is touching the skin, adding conductive gel to the electrodes if necessary.

**Step 4:** Treatment

Once the device is adjusted, you can move around easily.

**Auto-lock:** If you no longer press on the buttons, after 30 seconds the stimulation level increase function will lock to prevent an accidental increase. To increase the intensity, first unlock the device by pressing on the (-) button, then press on the (+) button to increase the intensity. Then let the programme run until it is complete.

The programme lasts 20 minutes. After 20 minutes, the device switches off automatically.

**Step 5:** Stop the device

The device can be switched off at any time by pressing and holding down the on/off button.

Do not remove the device from your skin during the treatment. Always turn off the device before removing it.

**Step 6:** Clean the device then put it away.

### b. Cleaning

The device must be switched off and unplugged for maintenance and cleaning.

Cleaning the electro-stimulator: Clean the device with a soft, smooth, dry cloth.

Cleaning the electrodes: After use, it is advisable to wipe the electrodes with a damp cloth to remove any remaining gel.


We recommend that you clean the electrodes at least once a week when used daily. Remove the electrodes before cleaning them. Clean the electrodes separately from the device, with clean, lukewarm water. Make sure they are clean and dry before repositioning them on the device.

Cleaning the strap: Handwash the strap separately. Avoid machine washing the strap.

### c. Charging the device

To charge the device, plug in the USB lead provided to the micro-USB port and the other end of the lead to a USB power socket (5V direct current) at least overnight (8 hours). A charge indicator will light up (next to the (+) button) and switch off when the battery is fully charged.

We advise you to fully charge the device at least once every 3 months to maximise battery life.

## 7. STORAGE CONDITIONS

Store the electro-stimulator, the electrodes and the manual in appropriate packaging in a dry place at room temperature (between -10°C and +40°C).

## 8. DISPOSAL OF THE DEVICE PARTS

Dispose of the device in accordance with Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Contact your local waste disposal authorities if in doubt. Electrodes and straps: these parts are not hazardous waste and can be disposed of with household waste.

## 9. WARRANTY

The device is guaranteed for 2 years from the date of receipt. The warranty does not cover the electrodes or any breakage or failure when the device is not used as instructed in the manual. Questions or complaints: contact the manufacturer by email at service.client@stimuli-technology.com or by post at the following address: Stimuli Technology, 20 bis rue Barthélémy Danjou, 92100 Boulogne-Billancourt, France, or seek advice from a healthcare professional.

## 10. TROUBLESHOOTING

| Situation  | Possible causes                                      | Solutions  |
|--|--|--|
| The treatment seems to be different or more unpleasant than before | The electrodes are not positioned in the right place | Switch off the device and move the electrodes slightly   |
|  | Intensity is too high or low                         | Change the intensity using the intensity level buttons   |
|  | The electrodes are not directly touching the skin    | Switch off the device and remove the device from the skin. Add a drop of conductive gel to the electrodes and reposition the device so that it is in contact with the skin |

| Situation  | Possible causes   | Solutions   |
|--|---|---|
| The treatment seems to be different or more unpleasant         | The electrodes are too dirty  | Clean the electrodes (section 8.b). If the problem occurs again, replace the electrodes |
|  | The electrodes are worn out   | Remove the electrodes and replace them with new ones (section 8.a)                      |
| The device turns off or stops the stimulation during treatment | The device turns off automatically if there is no activity or stimulation after 2 minutes | This is normal: switch the device back on   |
|  | The device turns off automatically after 20 minutes of treatment                          | This is normal: the treatment is finished   |
|  | The battery is too low  | Charge the device. You can use the device once it is charged                            |
|  | The device is no longer in contact with the skin or the device has come off               | Adjust the device again   |
| The device will not turn on                                    | The battery is too low  | Charge the battery  |
|  | The device is defective   | Contact customer service  |

## 12. ELECTROMAGNETIC SENSITIVITY

This device complies with the requirements of standard EN 60601-1-2 describing electromagnetic compatibility (EMC) conditions for medical devices. Tensi+ requires EMC precautions. Tensi+ must be installed and used in accordance with the EMC recommendations below.

Compliance with EMC standards does not mean that a device is totally immune. Tensi+ can be affected by portable or mobile RF communication equipment.

Tensi+ should not be used near to or stacked with other devices. If this is not possible, monitor the device to check that it operates normally under these conditions.

Risk of interference: the use of accessories and leads other than those specified, with the exception of leads and accessories sold by the manufacturer as spare parts for internal components, may increase the device's emission levels or decrease its immunity levels.

See electromagnetic emissions and immunity tables on the following pages.

The patient and/or user must notify the manufacturer, Stimuli Technology, of any serious incident involving the product at: *service.client@stimuli-technology.com* or Stimuli Technology, 20 bis rue Barthélémy Danjou, 92100 Boulogne-Billancourt, France, and notify the competent local authority.

| Guidance and manufacturer's declaration - electromagnetic emissions  |            |   |
|--|------------|---|
| Tensi+ is intended for use in the electromagnetic environment specified below. The patient or user should ensure that it is used in such an environment. |            |   |
| Emissions test   | Compliance | Electromagnetic environment - Comments                      |
| Electromagnetic radiation disturbance (Radiated emissions) (CISPR 11)  | Group 1    | The medical device uses RF energy for internal functioning. |
| Interference terminal voltage (conducted emissions) (CISPR 11)   | Class B    |   |

| Guidance and manufacturer's declaration - electromagnetic immunity   |                                     |                                     |   |
|--|-------------------------------------|-------------------------------------|---|
| Tensi+ is intended for use in the electromagnetic environment specified below. The patient or user should ensure that it is used in such an environment. |                                     |                                     |   |
| Immunity test  | IEC 60601 test level                | Level of compliance                 | Electromagnetic environment - Guidance                            |
| Electrostatic discharge (ESD) (IEC 61000-4-2)  | ± 8 kV on contact<br>± 15 kV in air | ± 8 kV on contact<br>± 15 kV in air | Home healthcare environment and professional healthcare facility. |
| Industrial frequency magnetic field (IEC 61000-4-8)  | 30 A/m                              | 30 A/m                              | Home healthcare environment and professional healthcare facility. |

| <b>Guidance and manufacturer's declaration -<br/>electromagnetic immunity, radiofrequencies</b>   |   |   |   |
|---|---|---|---|
| Tensi+ is intended for use in the electromagnetic environment specified below. The patient or user should ensure that it is used in such an environment.  |   |   |   |
| WARNING: Portable RF communication equipment (including peripherals such as antenna leads and external antennas) should not be used within 30 cm (12 inches) of any part of the TESTED DEVICE, including leads specified by the manufacturer. The performance of these devices could be impaired by such use. |   |   |   |
| <b>Immunity test</b>  | <b>IEC 60601 test level</b>   | <b>Compliance level</b>   | <b>Electromagnetic environment -<br/>Guidance</b>                 |
| Radiated radiofrequency electromagnetic fields (IEC 61000-4-3)  | 3 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz   | 3 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz   | Home healthcare environment                                       |
|   | 10 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz  | 10 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz  | Professional healthcare facility                                  |
| Nearby fields emitted by RF wireless communication devices (IEC 61000-4-3 provisional method)   | 9 V/m<br>710 MHz, 745 MHz,<br>780 MHz, 5240 MHz,<br>5550 MHz, 5785 MHz,                       | 9 V/m<br>710 MHz, 745 MHz,<br>780 MHz, 5240 MHz,<br>5550 MHz, 5785 MHz,                       | Home healthcare environment and professional healthcare facility. |
|   | 27 V/m<br>385 MHz   | 27 V/m<br>385 MHz   |   |
|   | 28 V/m<br>450 MHz, 810 MHz,<br>870 MHz, 930 MHz,<br>1720 MHz, 1845 MHz,<br>1970 MHz, 2450 MHz | 28 V/m<br>450 MHz, 810 MHz,<br>870 MHz, 930 MHz,<br>1720 MHz, 1845 MHz,<br>1970 MHz, 2450 MHz |   |

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## EN symbols glossary

|  |   |
|--|---|
|  | This manual contains information related to you safety.   |
|  | CE marking. This device meets the regulatory requirements of European regulation 2017/745 on medical devices.   |
|  | Device with type BF protection against electric shocks.   |
|  | Device serial number <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> |
|  | Part number   |
|  | Protection against foreign objects >12.5 mm. Water resistance: protection against vertically falling water drops (when tilted up to 15°).   |
|  | Name and address of manufacturer  |
|  | The device, its accessories and its packaging must be recycled appropriately when no longer used. Please follow local rules and regulations.  |
|  | This product is a medical device.   |
|  | Read the user manual carefully before using the device.   |
|  | Do not use the product if either the product or packaging is damaged.   |
|  | Storage conditions: keep away from moisture and store at between -10°C and 40°C   |
|  | Storage conditions: relative humidity 15 to 90%, Pressure 700 to 1060 hPa   |
|  | Not to be used with a pacemaker / cardiac electrical stimulation device   |

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