

上海大胜卫生用品制造有限公司

Shanghai Dasheng Health Products Manufacture Co., Ltd.

228 shihui Road, Songjiang, Shanghai 201613, China

Technical Datasheet

Model No.: FFP2 NRD

DTC3B, DTC3M, DTC3W, DTC3X, Small A-F, 2210, 2220, DAC4C

Ultrasonic welding technology, glue-free and odorless.

Arc-ends design of nose clip, avoid scratching the face, safer to wear

Material: PP Filtration Fabric (Melt Blown)

With earloops

With Nose Clip

Color: White

Bacterial Filtration Effectiveness $\geq 94\%$

Breath resistance: $(\Delta P) < 5$

Pressure difference: $(\Delta P) < 49.0$

