


Manual técnico

Vainas de succión- evacuación Clearpetra.

A) Vaina Acceso ureteral.


Dimensión de la vaina.

Tamaño	I.D. (+/- 0.3mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
10fr	3.3	4.0	13/26/36/40/46/55
11fr	3.7	4.3	
12fr	4.0	4.7	
13fr	4.3	5.0	
14fr	4.7	5.3	
Las vainas de las longitudes 13/26/36 tienen una salida presión. No así las medidas 40/46/55cm.			



Dimensión del obturador.


Tamaño	I.D. (mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
10fr	>/- 1.0	3.3	20/33/43/47/53/62
11fr		3.7	
12fr		4.0	
13fr		4.3	
14fr		4.7	



B) Vaina nefrostomía.


Dimensión de la vaina.

Tamaño	I.D. (+/- 0.3mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
10fr	3.3	4.0	13/17/21
12fr	4.0	4.7	
14fr	4.7	5.3	
16fr	5.3	6.0	
18fr	6.0	6.7	
20fr	6.7	7.3	
22fr	7.3	8.0	



Dimensión del obturador.


Tamaño	I.D. (mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
10fr	>/- 1.0	3.3	20/24/28
12fr		4.0	
14fr		4.7	
16fr		5.3	
18fr		6.0	
20fr		6.7	
22fr		7.3	



c) Vaina cistoscopio.

Dimensión de la vaina.

Tamaño	I.D. (+/- 0.3mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
18fr	6.0	6.7	21/24
20fr	6.7	7.3	
22fr	7.3	8.0	



Dimensión del obturador.

Tamaño	I.D. (mm)	O.D. (+/- 0.3mm)	L (+/- 2.5cm)
18fr	>/- 1.0	6.0	27/30
20fr		6.7	
22fr		7.3	

